



ANALYSIS INSTRUMENTS

WHERE SERVICE AND QUALITY IS OUR FIRST NATURE

Hand Crimper, Decrimper (MAKE : ANALYSIS)

Quality

- The vial crimper, mouth and handle part all is made of the stainless steel material, it is easy to clean and rustiness.
- Every tool made by us is customized, made one by one, not by production line.
- We has the ability to design and manufacture vial crimpers to meet special or unique requirements, so the cap sealing on bottle will be perfect, 100% matched.

Product Description

- Achieve secure, repeatable crimps or remove seals with minimized risk of damage.
- Lightweight, precision crimpers attach seals to chromatography vials. Available for standard aluminum or flip-off seals. Factory set to properly crimp standard-thickness seals. Adjustable handle setting screw allows proper crimps with septa of varying thickness.
- Forged pliers feature a fitted gripping surface designed to allow removal of the listed sizes of seals. Simply grip pliers around top of seal, squeeze and pry upward to remove.



Features

- Hand Crimper is widely used in laboratory Headspace vials crimping caps.
- The vial crimper is ideal for low volume crimping operations and sample preparation operations.
- The capping is safe, credible, tight, slick, good sealing.
- The vial crimper is ideal for low volume crimping operations and sample preparation operations.
- The capping is safe, credible, tight, slick, good sealing.
- Crimping machine reduced housing pumps offers a lower cost alternative to screw-on pumps, yet provides an upscale look to your finished package.
- Crimping machine is safe and practical, there is no requirement for electricity or compressed air, trouble-free operation.
- Crimping machine head is machined to precise limits, more accuracy.





ANALYSIS INSTRUMENTS

WHERE SERVICE AND QUALITY IS OUR FIRST NATURE

How to use

This machine is widely used in perfume industry, where small and medium scale production is required. If there is more production required, and if it is more human labour available, then also this machine is suitable.

- Put the spray on bottle, put the crimping head above perfume pump and hold the crimping hand tight.
- The tool covers the spray and a pulling action of the tool is occurring during pressing the handle. This causes the perfect cap sealing to the bottle.

Technical Parameter

- There are 2 parts crimping head and crimping hand.
- The hand crimper is designed and built with an adjustable pressure block in the crimper head that allows the operator to adjust the overall length of stroke that the cone moves when the seal is crimped.

MOC

- Stainless Steel
- MS Powder Coated

Guarantees

- Life-long repair.
- Quality problem (except for man-made damage), return and change, all is free
- Our customer service and technical department is available to answer your questions about our product or need help to decide which one is the right one for you.

Available Sizes listed as below

11 mm Crimper & Decrimper
13 mm Crimper & Decrimper
20 mm Crimper & Decrimper
32 mm Crimper & Decrimper





ANALYSIS INSTRUMENTS

WHERE SERVICE AND QUALITY IS OUR FIRST NATURE

Product Code:

Description	11 mm Crimper	Description	11 mm Decrimper
MOC	Stainless Steel/MS Powder Coated	MOC	Stainless Steel/MS Powder Coated
Part No.	AICRSS11/AICRMS11	Part No.	AIDCRSS11/AIDRMS11

Description	13 mm Crimper	Description	13 mm Decrimper
MOC	Stainless Steel/MS Powder Coated	MOC	Stainless Steel/MS Powder Coated
Part No.	AICRSS13/AICRMS13	Part No.	AIDCRSS13/AIDRMS13

Description	20 mm Crimper	Description	20 mm Decrimper
MOC	Stainless Steel/MS Powder Coated	MOC	Stainless Steel/MS Powder Coated
Part No.	AICRSS20/AICRMS20	Part No.	AIDCRSS20/AIDRMS20

Description	32 mm Crimper	Description	32 mm Decrimper
MOC	Stainless Steel/MS Powder Coated	MOC	Stainless Steel/MS Powder Coated
Part No.	AICRSS32/AICRMS32	Part No.	AIDCRSS32/AIDRMS32



AT info.analysisinstruments@gmail.com

ADD: SHOP NO. 10, SAI DEEP PLAZA,
NEAR D-MART, YASHWANT VIVA
TOWNSHIP, VASAI-VIRAR-401209
TEL : +91-9867226993



HEAD OFFICE:- SHOP NO.10.SAI DEEP PLAZA.YASHWANT VIVA TOWNSHIP.VASAI-VIRAR-401209
E MAIL:- info.analysisinstruments@gmail.com CONTACT NO: +91-9867226993